

Liquid-Liquid extraction pilot unit - rotative disc technology



Experimental capabilities

- Study of a liquid-liquid extraction pilot
- Study of rotating disk technology
- Study of the efficacy of the extraction column
- Study of mass balance
- Calculation of the number of theoretical stages

DIDATEC– Zone d'activité du parc – 42490 FRAISSES- FRANCE
Tél. +33(0)4.77.10.10.10 – Fax+33(0)4.77.61.56.49 – www.didatec-technologie.com
email : service_commercial@didatec-technologie.com

Reproduction interdite / copy prohibited– Copyright DIDATEC janv.-17- page 1

Dans le cadre de l'amélioration permanente de nos produits, ce descriptif technique est susceptible d'être modifié sans préavis
As part of the continuous improvement of our products, this technical specification may be modified without previous notifying

Operating principle

The GPC E11 bench allows the study of a liquid-liquid extraction system with rotating discs.

A gear pump provides supply of the mixture through the top of the column. A second gear pump ensures the solvent supply from the bottom of the column. The mixture and the solvent therefore will circulate at counter-current, which will allow the solvent to extract the acetic acid from the water by shearing of the discs on the liquid.



Students will study the liquid-liquid extraction system.

The robust design of this device makes it suitable for use in schools.

The equipment is set up on an Anodized aluminium frame on casters wheels. This gives it great strength and a flexibility of integration into your premises.

The manufacture of this equipment complies with the European standard for machinery manufacturing.

Illustrations



Technical details

1. **Agitator**
- Composed of a variable speed geared motor unit
2. **Extraction leg**
- Material: borosilicate glass
- Allows the interface level control
3. **Extraction column**
- Material: borosilicate glass
- DN : 50 mm
- Height: 2000 mm with 25 stages
4. **Digital flowmeters**
- Scale: 0-25 L/min
5. **Sample taking valve**
- Extract and raffinate
6. **Two supply bins**
- Material: polyethylene
- Volume : 35 L
7. **Electrical box**
8. **Two recovery bins**
- Material: polyethylene
- Volume : 35 L
9. **Glass decanter**
- Volume : 3 L

Two gear pumps
- With controllable variable speed drives

Temperature probes Pt100

- T1 : Base
- T2 : 1/3 lower of the column
- T3 : 1/3 higher of the column
- T4 : Column head
- T5 : Supply of the mixture
- T6 : Solvent supply

System of acquisition and computerized control

Control software and data acquisition

GPCE11



Services required

- Electrical supply : 230 Vac – 50 Hz – 20 A
- Electrical network : 1 phase) + Neutral +
- Dimensions: (L x W x H mm): 1800 x 760 x 2800
- weight (Kg): 200

Note : if the equipment installation is operated by our staff, all supplies and exhaust connections required must stand at less than 2m from the machine


Documentation

- User's manual
- Pedagogical manual
- Technical documentation of the components
- Lab exercises
- Software
- Certificate of conformity CE

Monitoring: Parameter setting, Plot of curve, control

The bench is also equipped as standard with a monitoring and configuration software. The connection towards the PC is made via a standard USB port. The software is divided into three parts:


SYNOPTIC:



We find in this window the block diagram of the machine with the location of various measures of process and their values.

GRAPH:

We find in this graph window, the possibility of drawing the measurement curves as a function of the time by selecting the desired quantities.



DIDATEC– Zone d'activité du parc – 42490 FRAISSES- FRANCE
Tél. +33(0)4.77.10.10.10 – Fax+33(0)4.77.61.56.49 – www.didatec-technologie.com
email : service_commercial@didatec-technologie.com

Reproduction interdite / copy prohibited– Copyright DIDATEC janv.-17- page 3

Dans le cadre de l'amélioration permanente de nos produits, ce descriptif technique est susceptible d'être modifié sans préavis
As part of the continuous improvement of our products, this technical specification may be modified without previous notifying

Illustrations non contractuelles / Illustrations not contractual

version : FT-GPCE11-STD-A

GPCE11



CONTROL-PARAMETER:

We find in this tab, the possibility to control via the supervision the powering and adjusting the speed of the pump 1, pump 2 and agitator

